Návod na použitie

MESI mTABLET

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🆕 MESI	16:32	24th of July 2016 🜩	92% 🧰	6 MESI 16:32	24th of July 2016 💎 92% 🧰
Choose Wi-Fi			_	Register I am a new user.	
Hospital Patient - Secured ✓ MESI Wi-Fi Staff Public - Guest Dr. Smith's iPhone		8 8 8 8	* * * *	E-mail Enter e-mail address Phone number Enter your phone number Country Germany	
BACK		+Add net		Register Description Login I have mRecords account. Item Records account.	

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Locace agreement for Mesi cloud Section 2 (Mesi Cloud) These terms of use (these Terms') apply to visitors of MESI's ('we' 'us' or 'our') webaits meticia, www.meiamedical.com, and mecords meiaid, as well as registered users of our Mesives). Please review these Terms carefully. Please also review on Princey Policy whence years enview these Terms carefully. Please also review on Princey Policy sections a device of the section of the sec	ESI	1	6:32	24th of July 2016 🛛 😤	92%
Even terms of use (these "Terms") apply to visitors of MESI's ('we' 'us' or 'our') websites metais, www.meannedical.com, and intercods meiais, as well as registered users of our MESI's - and mRecords - information and technology services and mobile applications (our 'Services') Please relevent these Terms and talk and protect you privacy Policy, which explains how our Services they are personal data and protect you privacy when you are consistent and your personal data and protect you privacy. Policy when you are consistent and your personal data and protect you privacy Policy when you are consistent and your personal data and protect you privacy. Policy when you are consistent and your personal data and protect you privacy Policy when you are consistent and your personal data and protect you privacy. Policy when you are consistent and your personal data and protect your privacy Policy when you are consistent and your personal use and protect your privacy. Policy services are allocide to the policic Services are appreciable to the policic Services on long as you comply with trees Terms. This personal use includes linking to information appearing in the Public Services. Proved and the protect is an endipole, agent or regress yours are and the public Services are on the protect to private and the public Services are on a privacy servet.					
Exercise agreement for Mesi cloud Exercise terms of use (these "Terms") apply to visitors of MESI's (we" vis" or 'our") website mesia, www.meemendeal.com, and mercercist, mesia, well as registered users of our MESI®- and mRecords®-information and technology services and mobile applications (our Sarvices). Please review these Terms and etchnology services and mobile applications (our Sarvices). Please review these Terms and etchnology services and mobile applications (our Sarvices). Please review these Terms and etchnology services and mobile applications (our Sarvices). Please review these Terms and etchnology services and mobile applications (our Sarvices). Please review these Terms and our Privacy Policy, which explains how car's devices travelup to previous the terms and our Privacy Policy. I. Accessing our Services. I. Dublic Services Thublic Services Thuble Services. Thuble Services available without registering or obtaining a password. We call these Thuble Services. You may make personal, non-commercial use of the Public Services so long as you comply with these Terms. This personal use includes linking to information appearing in the Puble Services. This personal use ado includes temporary caching by you browsers or a proy serve. If you perste a free, public search engine that has made a public commitment to adhering to the robots to protocol. Longether with our Stemapa and Crawdeley directive (collective), the Vistor of the provision or display of aid hyperlinks (including so called 'deep linki') to the Public Services or on gas (b) you do so is an orper provision or adiative of you is and cached pages, and (i) if our user again it adiatived by puit and vision of instructions public services notify and (i) in public Services, so long as (b) you do so is an orper provision and cached pages, and (ii) if our user again it adiative data you is an orper term for there is a service so and (ii) if our u					
These terms of use (these "Terms") apply to visitors of MESI's ("we" vas" or "our") websites mei si, www.meimedical.com, and mrecords mei si, as well as registered users of our MESI's and mrecords- information and technology reviews and mobile septilations (our "brivices"). Please review these Terms carefully. Please also review our Privacy Policy, which explains how our Services, the areases in our services and mobile septilations (our Services). Accessing our Services were aprecedent to these Terms and our Privacy Policy, which explains how our Services and mobile services and mobile septilations (our Services). Baces visces the accessing our Services and mobile services are provided with the relations to information appearing in the Public Services. Provided you do not represent visours of as an employee, agent or representative of MESI, Ltd. Personal use also includes Imising to information appearing the Public Services and our Services, (i) you do not directly or carbing by your browser or a provide of MESI, Ltd. Personal use also includes temporary caching by your browser or a protocol), you may cardii do and publish hyperitiks (including so called 'deep links') to the Protocol motile with the Protocol instructions public and publish hyperitiks (including so called 'deep links') to the Protocol on services, (i) you do not directly or informations page. and find part user-agent is disallowed by us in our Protocol instructions or we otherwise notify you in writing you shall cases all carving of our Services. We call such a free, public search engine statifying the requirements of clauses (i) through (ii), a "Pu	E Licence a	greement for Mesi			
These terms of use (these "Terms") apply to visitors of MESI's (twe" vis" or 'our") websites meisia, www.meismediaal.com, and mreccords meisia, as well as registered users of our MESI's and mRecords-information and technology services and mobile applications (our "Sirvices"). Please review these Terms and rully. Please also review our Privacy Policy which explains how our Services trans carefully. Please also review our Privacy Policy which explains how our Services and mobile applications (our Sirvices). Please review these Terms and early privacy. These methods explains how our Services are your personal data and protect your privacy when you use our Services. By accessing our Services, you are agreeing to these Terms and our Privacy Policy. The Public Services are also includes the privacy and the public Services are applied by a privacy. These terms are privacy and the privacy and the public Services are one of MESI. Ltd. Personal use also includes temporary caching by your browser or a privacy service. They caprate free, public Services are one one one please and the Public Services are one of MESI. Ltd. Personal use also includes temporary caching by your browser or a privacy service. They caprate free, public Services one one one one tempore with the Protocol. Together with our Stemaps and Craw-delay directive (collectively, the Protocol), you may craw. Index and public the proteins on devices (privations publiced in applicable locations on or Services. (b) you do not directly or indexetly receive returners ion in connection with the provision or display of alsi hypetinks and cached pages, and (ii) of your user agent is disallowed by usi in our Protocol instructions public service tengine and the public Services. Services Weall such a free, public search engine and the public Services. We call such a free, public search engine and write you shall cases all crawling of our Services. We call such a free, public search engine and write you shall cases all crawling of our Services. We call such a free, publi	_				
Interaction to the problem of the proportion of the problem of the	These terms of use (these "Terms") apply to y	isitors of MESI's (*)	we' 'us' or 'our') websites	
MESIA- and mRecordsB- information and technology services and mobile applications (our Sarvices). Please review these Terms carefully, Please also review one Privacy Policy, which explains how careful services transport personal data and protect you privacy when you are our Sarvices. By accessing our Services, you are agreeing to these Terms and our Privacy Policy. J. Accessing our Services . Burden Sarvices and allable without registering or obtaining a password. We call these "Public Services" You may make personal, non-commercial use of the Public Services so long as you comply with these Terms. This personal use includes linking to information appearing in the Public Services. If you not net represent yourself as an employee, agent or repre- sentative of MESI, Ltd. Personal use also includes temporary caching by your browser or a proxy server.	mesi.si, www.mesim	edical.com, and mrecord	s.mesi.si, as well a	is registered users of our	
Stavies I) Please review these Terms carefully Please also review our Privacy Policy, which explains how cur Services it year personal data and protect your privacy when you use our Services it ya accessing our Services, you are agreeing to these Terms and our Privacy Policy. 1. Accessing our Services 3.1 Public Services We make some Services available without registering or obtaining a password. We call these Public Services 'You may make personal use also includes linking to information appearing in the Public Services. Provided you do not represent yourself as an employee, agent or representative of MESL tot. Personal use also includes linking to information appearing in the Public Services, provided you do not represent yourself as an employee, agent or representative of MESL tot. Personal use also includes linking to information appearing to the robots to protocol, together with our Stemage and Craw-Gelg vincetwy, the 'Potocol', you may crawl, index and publich hyperlinki (including so called 'deep linking') to the 'Potocol', you may crawl, index and publich hyperlinkis (including so called 'deep linking') to the 'Potocol', you may crawl, index and publich hyperlinkis (including so called 'deep linking') to the 'Potocol', you may crawl, index and publich hyperlinkis (including so called 'deep linking') to the 'Potocol', you shall ceaser all crawling of our Services. We call such a free, public search engine statisty ing the requirements of clauses () through (ii), a "Public Service tengine' and the Potoloc instructions publicies each engine' antistry put shall ceaser all crawling of our Services. We call such a free, public search engine' and the Public search engine' and th	MESI®- and mRecord	ds®- information and tec	hnology services a	and mobile applications (c	ur
explains how out Bevices treat your personal data and protect your privacy when you use our Services. By accessing our Services, you are agreeing to these Terms and our Privacy Policy. 1. Accessing our Services 1.1 Public Services We make some Services available without registering or obtaining a password. We call these "Public Services available without registering or obtaining a password. We call these "Public Services available without registering or obtaining a password. We call these "Public Services available without registering or obtaining a password. We call these "Public Services provided you do not represent yourself as an employee, agent or repre- sentative of MESI, Ltd. Personal use also includes temporary caching by your browner or a proy serve: Hyou operate free, public search engine that has made a public commitment to adhering to the robots tht protocol. together with our Sitemaps and Crawi-delay directive (collectively, the "Pholic Services long as (i) you do not directly or indirectly receive remurses- tion in connection with the provision or display of aid hypeinks inducides pages, and (ii) of your user-agent is disallowed by us in our Protocol instructions or we deriveise notify you writing you shall cease all crawing of our Services. We call such a free, public search engine satisfying the requirements of clauses (i) through (ii), a "Public Search Engine" and the Public BACK AGREE	"Services"). Please re	view these Terms carefu	lly. Please also rev	iew our Privacy Policy, wh	ich
Services By accessing our Services, you are agreeing to these Terms and our Phixacy Policy: A. Accessing our Services A. Dublic Services TJ Public Services We make some Services available without registering or obtaining a password. We call these Trublic Services 'You may make personal, non commercial use of the Public Services so long as you comply with these Terms. This personal use includes linking to information appending in the Public Services.' If you operate a free, public search engine that has made a public commitment to adhering to the robots to provise and on the resenser yourself as an englose, agent or rapper sentative of MESI, Ltd. Personal use also includes temporary caching by your browser or a proxy serve. If you operate a free, public search engine that has made a public commitment to adhering to the robots to protocol. together with our Stemapa and Crawfeley directive (collective), the 'Photocol', you may craw, index and public Microticons or we otherwise notify you use any one is one compliance with the Protocol instructions public Services one age (9) you do not directly or indirectly receive remuners- tion in connection with the provision of adaptor of aid hypelinks and cached pages, and (ii) if your user-agent is disaltword by us in our Protocol instructions or we otherwise notify you on the requirements of clauses (i) through (ii), a "Public Search Engine", and the Public	explains how our Ser	vices treat your personal	data and protect y	your privacy when you use	our
1. Accessing our Services House Some Services wallable without registering or obtaining a password. We call these "Public Services vious on the personal, non-commercial use of the Public Services so long as you comply with these Terms. This personal use includes linking to information appearing in the Public Services, provided you do not regressent yourself as an employee, agent or regressentative of MEX. It. A Personal use also includes temporary caching by your browser or a proxy serve. If you operate a free, public search engine that has made a public commitment to adhering to the robotic stop protocol, together with our Stemapa and Craw delay directive (collectively, the "Protocol"), your grawn, index and publich hyperitiking includings a called deep linking to the Public Services, so long as (i) you do so in compliance with the Protocol instructions published in appliciable locations on our Services, (ii) you do not directly or indirectly receive remuners to in connection with the provision or display of and thyperlinks and cached pages, and (iii) provide search engine startisty of MEX. BACK Agree	Services. By accessir	ng our Services, you are a	igreeing to these T	erms and our Privacy Poli	CY.
1.1 Public Services We make some Services available without registering or obtaining a password. We call these public Services are large some the Public Services so long as you comply with these Terms. This personal use includes linking the individual services are provided you do not represent yourself as an employee, agent or rapresentative of MES. Ut Personal use also includes temporary caching by your browser or a proxy serve. If you operate a free, public search engine that has made a public commitment to adhering to the robotic sup rotocol, together with our Stempag and Craw delay directive (collectively, the "Potocol"), your your your, index and publich hypetifiks (including so called be pass, and (ing) is not compliance with the Protocol instructions publicable locations on our Services. (i) you do not directly or indirectly receive remainers to in connection with the provision of splayd or add hypetifiks and cabet pages, and (ing) you do search engine that we applie link and cabet pages. The provision is uncertained to the public services we call such as free, public search engine that has made a public commitment to adhering to the protocol instructions on we otherwise notify you in writing you shall cease all crawfing of our Services. We call such a free, public search engine, and the Public service search engine, and the Public search e	1. Accessing	our Services			
We make some Services available without registering or obtaining a password. We call these "Public Services are not an expressional non-commercial use of the Public Services so long as you comply with these Terms. This personal use includes linking to information appending in the Public Services, provided you do not represent yourself as an employee, agent or repre- sentative of MESI, Ltd. Personal use also includes temporary caching by your browser or a provy server. If you operate a free, public search engine that has made a public commitment to adhering to the robots to portice. They prove the protocol, together with our Sitemapa and Crawf-delay directive (collectively, the "Photocol"), you may craw, index and public hyperlinks (including so-called 'deep links') to the Photos Sories as (b) you do not directly or inderectly receive remurses- tion in connection with the provision or display of said hyperlinks and cached pages, and (iii) if your user-agent is disallowed by usi in our Protocol instructions public services pages, and (iii) if your user-agent is disallowed by usi in our Protocol instructions or we otherwise notify you in writing you shall caese all crawling of our Services. We call such a free, public search engine satisfying the requirements of clauses (i) through (ii), a "Public Search Engine" and the Public BACK AGREE	1.1 Public Servi	ices			
Trubic Services 'You may make personal, non-commercial use of the Public Services so long as you comply with these Terms. This personal use includes linking to information appending in the Public Services, pervided you on to represent yourself as a menipower, agent or rapper sentative of MESI, Ltd. Personal use also includes temporary caching by your browser or a proxy server. If you operate a free, public search engine that has made a public commitment to adhering to the robots to protocol. together with our Stemapa and Crawid eley directive (collective), the 'Pototocol', you may crawi, index and publich hyperlinks (including so called 'deep links') to the Photosol', you may crawi, index and publich hyperlinks (including so called 'deep links') to the Photosol', you one craw on expression of display of said hyperlinks and cached pages, and (ii) if your user-agent is disaltworth by un on Photocol Instructions publicher ensures tori in connection with the provision of display of said hyperlinks and cached pages. and (ii) if your user-agent is disaltworth by un on Photocol Instructions or we otherwise notify you writing you shall cease all crawling of our Services. We call such a free, public search engine satisfying the requirements of clauses (i) through (ii), a "Public Search Engine", and the Public BACK AGREE	We make some Servi	ices available without req	istering or obtaining	ng a password. We call the	ese
as you comply with these Terms. This personal use includes linking to information appearing in the Public Services, provided you do not represent yourself as an employee, agent or repre- sentative of MES, Lit Dersonal use also includes temporary caching by your browser or a proxy server. If you operate a free, public search engine that has made a public commitment to adhering to the robots to protocol, together with our Stempag and Craw delay directive (collectively, the Photocol), you may craw, index and uplicits hyperlinks (includings or called deep links) to the Public Services, so long as (i) you do so in compliance with the Protocol instructions published in applicable locations on our Services, (ii) you do not directly or indirectly receive remanes- tion in connection with the provision of display of add hyperlinks and cached pages, and (iii) your user sagent is disallowed by us in our Protocol instructions or we otherwise notify you in writing, you shall cease all crawling of our Services. We call such a free, public search engine satisfying the requirements of clauses (i) through (ii), a "Public Search Engine" and the Public BACK Agree	"Public Services." You	u may make personal, no	n-commercial use	of the Public Services so I	ong
In the Public Services, provided you do not represent yourself as an employee, agent or repre- sentative of MSSL Ltd. Personal use also includes temporary caching by your browser or a proxy server. If you operate a free, public search engine that has made a public commitment to adhering to the robots to protocol, together with our Stemaps and Craw-Gelly directive (collectively, the "Protocol"), you may crawl, index and publich hyperlinks (including so called "deep links") to the Public Services, so long as (i) you do so in compliance with the Protocol instructions publiched in applicable locations on our Services, (i) you do not directly or infractly receive remuner- tion in connection with the provision of digaldy of aid hyperlinks and cached pages, and (iii) your user-agent is disallowed by us in our Protocol instructions or we otherwise notify you in writing you shall cease all crawling of our Services. We call such a free, public search engine satisfying the requirements of clauses (i) through (ii), a "Public Search Engine" and the Public BACK AGREE	as you comply with t	hese Terms. This person	al use includes link	king to information appear	ing
sentative of MESI, Ltd. Personal use also includes temporary caching by your browser or a proxy server. If you operate a free, public search engine that has made a public commitment to adhering to the robots to protocol. together with our Stemapa and Crawi-delay directive (collectively, the "Protocol"), you may crawi, index and publich hyperlinks (including so called 'deep links') to the Photos Services solong as (b) you do not directly or indirectly receive remunes- tion in connection with the provision or display of said hyperlinks and cached pages, and (iii) if your user-agent is diallowed by usi no 'Protocol instructions or wite therwise notify you writing you shall cease all crawing of our Services. We call such a free, public search engine satisfying the requirements of clauses (i) through (ii), a "Public Search Engine" and the Public BACK AGREE	in the Public Services	s, provided you do not rep	resent yourself as	an employee, agent or rep	pre-
If you operate a free, public search engine that has made a public commitment to adhering to the robots th protocol. Logether with our Sitemaps and Crawi-delay directive (collective); the "Potocol", you may crawi, index and publish hyperinks (including so-called 'deep links') to the Public Services, so long as () you do so in compliance with the Protocol instructions published in applicable locations on our Services, (ii) you do not directly receive remunes- tion in connection with the provision or display of said hyperlinks and cached page, and (iii) fragment satisfying the requirements of clauses (i) through (ii), a "Public Search Engine," and the Public BACK AGREE	sentative of MESI, Ltd	d. Personal use also inclu	ides temporary ca	ching by your browser or a	3
If you operate a free, public search engine that has made a public commitment to adhering to the robots to protocol, together with our Stemaps and Craw-delay directive (collectively, the "Potocol"), you may crawl, index and uplicit-hyperinks (includings ocalitadir deep links") to the Public Services, so long as (i) you do so in compliance with the Protocol instructions published in applicable locations on our Services, (i) you do not directly or indirectly receive remuners ton in connection with the provision of display of aid hyperlinks and cacher dapages, and (iii) your user-agent is disallowed by us in our Protocol instructions or we otherwise notify you in writing, you shall cease all crawfing of our Services. We call such a free, public search engine satisfying the requirements of clauses (i) through (ii), a "Public Search Engine" and the Public BACK AGREE	proxy server.				
the robots bit protocol, together with our Sitemaps and Crawi-delay directive (collectively, the Photocol), you may crawi-index and publish hyperfinits (including so-called 'deep links') to the Photol Services. So long as (b) you do is in compliance with the Protocol inductions published in applicable locations on our Services. (b) you do not directly or indirectly receive remunes- tion in connection with the provision or display of aaid hypefinits and cached pages, and (iii) if your user-agent is disallowed by us in our Protocol instructions or we otherwise notify you in writing you shall cease all crawling of our Services. We call such a free, public search engine satisfying the requirements of clauses (i) through (ii), a "Public Search Engine" and the Public BACK AGREE	If you operate a free,	public search engine that	t has made a publi	ic commitment to adherin	g to
Photocoff, you may craw Ladex and publish hyperlinks (including so called 'deep links') to the Public Services, so long as (i) you do so h compliance with the Protocol instructions published in applicable locations on our Services, (i) you do not directly or indirectly receive remanes- ton in connection with the provision of display of said hyperlinks and cached pages, and (iii) if your user-agent is disallowed by us in our Protocol instructions or we therewere notify you writing, you shall cesse all crawling of our Services. We call such a free, public search engine satisfying the requirements of clauses (i) through (iii), a "Public Search Engine", and the Public BACK AGREE	the robots bit protoco	ol, together with our Siter	maps and Crawl-de	elay directive (collectively,	the
Public Services, so long as () you do so in compliance with the Protocol instructions published in applicable locations on our Services, (ii) you do not directly receive remuners. tion in connection with the provision or display of said hypefinks and cached page, and (iii) fry your user agent is disallowed by us in our Protocol instructions or we otherwise notify you in writing, you shall ceesal al crawing or our Services. We call such as free, public search engine satisfying the requirements of clauses (i) through (ii), a "Public Search Engine," and the Public BACK AGREE	"Protocol"), you may	crawl, index and publish I	hyperlinks (includir	ng so-called "deep links") t	o the
In applicable locations on our Services, (i) you do not directly or indirectly receive remuner- tion in connection with the provision of odplayd of aid hypenlika sind cadebad pages, and (iii) if your user-agent is disallowed by us in our Protocol instructions or we otherwise notify you in writing you shall cease all crawing of our Services. We call such a free, public search engine satisfying the requirements of clauses (i) through (ii), a "Public Search Engine" and the Public BACK AGREE	Public Services, so lo	ung as (i) you do so in cor	npliance with the P	Protocol instructions publi	shed
too in connection with the provision of display of said hypefinitia and cached pages, and (iii) if your user-agent is disallowed by usin our Protocol instructions or we otherwise notify you in writing, you shall cesse all crawling of our Services. We call such a free, public search engine satisfying the requirements of clauses (i) through (iii), a "Public Search Engine," and the Public BACK AGREE	in applicable location	is on our Services, (ii) you	do not directly or	indirectly receive remune	ra-
boo lase span is a balaxies by due in our instructions of the outer the outer there interfy your in- writing you shall be seen all crawing on our Services. We due all such a free, public search Engine," and the Public back BACK AGREE	tion in connection will	th the provision or displa	y of said hyperlink:	s and cached pages, and ((iii) if
satisfying the requirements of clauses () through (ii), a "Public Search Engine," and the Public BACK AGREE	writing you shall ge	salowed by us in our Pro	vices. We call such	b a free public search eoo	ine
	satisfying the require	ments of clauses (i) thro	ugh (iii), a "Public \$	Search Engine," and the Pu	blic
	616.5		900		
		BACK	AG	REE	
	_				
		÷	0	0 5	
					n i
			0	(÷]

Registration	
*Name	Enter Name
*Sumame	Enter Surame
*Specialization	Choose
*Work sector	Choose
*Password	
*Repeat password	
This is also your password for mRecords	Show characters
*Un-lock code	
	Show characters
*Required fields	
BACK	NEXT
ല്പ	

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Sec. MESI	16:32 24th of July 2016 🍷 92% 🚥	● MESI 16:32 24th of July 2016 🖤 92% 🛲
Working group		Add users to working group
We have created a new w	orking group Smith Gregory, MD.	Add your coworkers to your working group to share your patient's data. Insert e-mail addresses of people you would like to invite.
Practice name	Enter practice name	You can manage users in Doctor > Settings > Working group.
Practice address		Enter e-mail addresses, separated by space
	Code City	
	Choose country 🗸	
Practice email		SKIP SEND
Practice website	URL address	
Practice phone		
Practice fax		
(ок	
	0	
	• • • •	

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Connect new module	Connect new module
mTablet has established connection with the right arm cuff module. To confirm it, press again the button on top of the module. If this is not the module you want to pair, tap BACK	Congratulations! You have successfully connected the right arm cuff module.
ВАСК	PAIR NEW FINISH

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Connect new module Trablet has established connection with the right arm ouff module. To confirm it, press again the button on top of the module. If this is not the module you want to pair, tap BACK.	<image/> <section-header><section-header><text><text></text></text></section-header></section-header>
ВАСК	PAIR NEW FINISH

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Decks Smith John MD Smith Gregory, GP Smith Gregory, GP Answer No patient selected Image: Answer Smith Gregory, GP Smith Gregory, GP Artications Artications Artications Artications Artications Smith Gregory, GP Smith Gregory, GP Artications Artications Artications Artications Smith Gregory, GP Smith Gregory, GP Smith Gregory, GP Artications Artications Artications Artications Smith Gregory, GP Smith Gregory, GP Smith Gregory, GP Artications Artications Artications Artications New patient Smith Gregory, GP Smith Grego	Image: Sinith John MD Image: Sinith John	S MESI	16:32	12.349 2014 🍄 213 82% 🗰	6 MESI	16:32	2400 of July 2016 💎 1929, 🚥	6 MESI	16:32	24th of July 2016 🍸 12% 🗰
Mexicanous Descurrence Water New patient Mexicanous Able ECG BP Patient data Oracle Office New patient Mexicanous Able ECG BP Bod Patient Oracle Office New patient Mexicanous Able ECG BP Bod Bernare Marrane <	ANDER No patient selected Image: And a	росток	Smith John MD Smith John's working group	\otimes	ростоя	Smith Gregory	y, GP	С остоя	Smith Gregor	y, GP
Arstantes Arstantes ECG BP Active Orarts Pedia Active Orarts Pedia Active Orarts Pedia Active Orarts Pedia BP Bool Pressore Pedia Break December 100 December 100 Active Pedia December 100 Break December 100 December 100 Break December 100 </td <td>APVECATION Archive Charits Poole ABI Mode Panare BP Bod Penare BD Bod Pen</td> <td>PATIENT</td> <td>No patient selec</td> <td>ted 🕀</td> <td>PATIENT</td> <td>New patient</td> <td>\otimes</td> <td>PATIENT</td> <td>New patient</td> <td>\otimes</td>	APVECATION Archive Charits Poole ABI Mode Panare BP Bod Penare BD Bod Pen	PATIENT	No patient selec	ted 🕀	PATIENT	New patient	\otimes	PATIENT	New patient	\otimes
ABI Acke buchti ECG Rectry Betrosendogen BP Buch hesser BUCh hesser BP Buch hesser <th< td=""><td>ABI Memory Note: ECC Memory ecosore/dogram BP Ded Pessare Partice data Concerner Manage Manage Mana</td><td></td><td></td><td></td><td>Archive</td><td>Charts</td><td>Profile</td><td>Archive</td><td>Charts</td><td>Profile</td></th<>	ABI Memory Note: ECC Memory ecosore/dogram BP Ded Pessare Partice data Concerner Manage Manage Mana				Archive	Charts	Profile	Archive	Charts	Profile
ABI ECG BP BP Thank the same t	ABI Add bachal Add bachal A				Patient data		⊘ ⊗	Patert o	lata	08
And bushel Percey Dec and any and any and any angle angle and angle a	And Buchal Petring Bood Pressure Bood Pressure Contact and the Boo	ABI	ECG	BP	Harre	*Sumame	*Race	"Name	*Sumame	*Race
Ansatz State of processing S	Alternative Autors Fundational Autors Fundational Fundational <td>Ankle Brachial Index</td> <td>Resting electrocardiogram</td> <td>Blood Pressure</td> <td>Enter name</td> <td>enter sumarne</td> <td>Caucasan 🗸</td> <td>Enter name</td> <td>Enter sumarne</td> <td>Caucasian 🗸</td>	Ankle Brachial Index	Resting electrocardiogram	Blood Pressure	Enter name	enter sumarne	Caucasan 🗸	Enter name	Enter sumarne	Caucasian 🗸
Heart No. Heart No. Address Prone number Email Address Prone number Email Direct Number Email Address Direct Street Number Email Address Color Ory Direct Direct Direct Direct Direct Direct	Heaver tess Heaver tess Heaver tess Address Phore number Ernal Address Phore number Dom Ofly Domal Address Domal Address Bornel Address Domal Address Color Ofly Domal Address Domal Address Domal Address Domal Address				Insurance number	DOWN/YYYY	Female: Male	Insurance number	DO/MM/1111	Female Male
Address Proce runber Email Address Proce runber Email Cover	Address Proce Auchor Ernal Dave Day Cone Day				Hequired fields			Pequired fields		
Strett Number Email Address Strett Number Email Address Code Ory Email Address	Dirett Number Dinal Address Dirett Number Diret Address Dirett Number Diret Address Dirett Number Dirett Address				Address	Phone number	E-mail	Address	Phone number	Email
	Covey Covey				Street Code City	Number	Email Address	Street Code City	Number	Email Address
Country Country					Country			Country		

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